


# VivaLAN™ – Network

Are you using different imaging systems in your clinic or practice? Are you tired of typing in the same patient information on each different device and of complicated backups?

Then the VivaLAN™ network solution is what you are looking for: imagine one system where all patient data is managed centrally and one software interface for clinical, dermoscopic and confocal imaging. All acquired images can be viewed on every computer connected to the network. A system that will help you to create a workflow that saves precious time and lets you work efficiently.



Local IT / Server  
Storage & BackUp



## ■ Advantages at a Glance

- Efficient workflow & data management
- Cost reduction through simplified data management
- Same software interface for clinical dermoscopic & confocal examination and documentation
- Optional connection to patient management systems (PACS, GDT and HL7)

## ■ Local IT Requirements

- Operating system      Microsoft Windows Server
- Database                Microsoft SQL Server
- Free storage space      2TB (seamless capacity expansion recommended)  
  ≈ 1 GB / patient
- Network connection    Gigabit Ethernet  
  *(A reliable network connection to the VivaScope System(s)  
  is required, no Wi-Fi)*

### VivaScope GmbH

Stahlgruberring 5  
81829 Munich, Germany  
Phone +49 (0) 89 / 420 96 - 280  
Email [info@vivascope.eu](mailto:info@vivascope.eu)

[www.vivascope.eu](http://www.vivascope.eu)

**VIVASCOPE**